

Taq DNA Polymerase 2x Premix with Dye



Catalog #	3249	3250
Package Size	500 reactions	1,000 reactions
Volume	12.5 ml	25 ml

Description

Intact Genomics (ig®) Taq DNA Polymerase 2x Premix with Dye is a ready to use premix which contains Taq DNA polymerase, dNTPs, MgCl₂ and stabilizers with optimized reaction buffer. It has been optimized for routine PCR applications. Taq is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity (1, 2) and a 5' flap endonuclease activity (3, 4). This product is supplied with the unique Intact Genomics 5x Magic Enhancer that enables efficient amplification of GC rich templates up to 84%.

Applications

- Routine PCR and RT-PCR
- Primer extension
- Colony PCR
- Genotyping
- Efficient for amplifying high GC content template DNA with Magic Enhancer.

Product Includes

- 1) Taq DNA Polymerase 2x Premix with Dye
- 2) 5x Magic Enhancer

Storage Temperature

-20 °C

1x Premix Composition

10 mM Tris-HCl pH 9.0
50 mM KCl
1.5 mM MgCl₂

0.2 mM dNTPs
5% Glycerol
0.08% Igepal CA 630
0.05% Tween-20
50 Units/ml Taq Polymerase.

Protocol

1. Thaw primer solutions, 5x Magic Enhancer (if required) and mix thoroughly before use.
2. Prepare a reaction mix according to the following table:

PCR Reaction Set Up:	
Template	1-50 ng
Forward Primer (3.2 µM)	1.0 µl
Reverse Primer (3.2 µM)	1.0 µl
5x Magic Enhancer (optional)	4.0 µl
Taq DNA Polymerase 2x Premix with Dye	10.0 µl
H ₂ O up to	20.0 µl

3. Mix the reaction mixture thoroughly.
4. Add template DNA to the individual PCR tubes containing the reaction mixture.
5. Program the thermal cycler according to the manufacturer's instructions. A typical PCR cycling program is outlined in the following table:

PCR Cycling Conditions			
Steps	Temp.	Time	Cycles
Initial Denaturation	95 °C	3 min	1
Denaturation	95 °C	30 sec	25-35
Annealing	55-60 °C	30 sec	
Extension	72 °C	1 min/kb	
Final Extension	72 °C	5-10 min	1
Hold	4 °C	∞	

6. Place the PCR tubes in the thermal cycler and start the cycling program.
7. Analyze 5 µl of PCR products by agarose gel electrophoresis.

Reference

1. Chien, A., Edgar, D.B. and Trela, J.M. (1976). J. Bact. 127, 1550-1557.
2. Lawyer, F.C. et al. (1993). PCR Methods and Appl. 2, 275-287.
3. Longley, M.J., Bennett, S.E. and Mosbaugh D.W. (1990). Nucleic Acids Res. 18, 7317-7322.
4. Lyamichev, V., Brow, M.A. and Dahlberg, J.E. (1993). Science. 260, 778-783.

Related Products

1. Hot Start Taq DNA Polymerase (Cat.# 3293)
2. Taq DNA Polymerase (Cat.# 3243)
3. i7® High-Fidelity DNA Polymerase 2X Master Mix (Cat.# 3257, 3259)
4. i7® Hot Start High-Fidelity DNA Polymerase 2X Master Mix (Cat.# 3284, 3286)

Technical Support

Intact Genomics is committed to supporting the worldwide scientific research community by supplying the highest quality reagents. Each new lot of our products is tested to ensure they meet the quality standards and specifications designated for the product. Please follow the instructions carefully and contact us if additional assistance is needed. We appreciate your business and your feedback regarding the performance of our products in your applications.

