



Glycerol-Free T4 UvsX DNA Recombinase

Manual

Catalog #	3562GF	3565GF	3567GF
Package Size	100µg	500µg	1000µg
Volume	20µl	100µl	200µl
Concentration	5µg/µl		



Important!

-20°C Storage Required

- * Immediately inspect packages
- * Freeze upon receipt



visit us online for more
products & custom services

Intact Genomics, Inc.

Table of Contents

Product Description.....	3
Protein Purity.....	3
Specifications.....	3
Components and Storage.....	3
Additional Buffer Information.....	3
Quality Control.....	4
Notes.....	4
Enzyme Concentration.....	4,5
Related Products.....	6
Ordering Information.....	6
References.....	6
Technical Support.....	7

Description:

Homologous recombination is important for the error-free repair of DNA double-strand breaks and for replication fork restart. Recombinases of the RecA/RAD51 family perform the central catalytic role in this process. UvsX recombinase is the RecA/Rad51 ortholog of bacteriophage T4. Intact Genomics (ig®) UvsX and other recombinases form presynaptic filaments on ssDNA that are activated to search for homology in dsDNA and to perform DNA strand exchange (1-3).

Protein Purity:

The physical purity of this enzyme is $\geq 98\%$ as assessed by SDS-PAGE with Coomassie® blue staining (Fig. 1).

Specifications:

Product Source: *E. coli* BL21 (DE3) strain expressing T4 UvsX gene

Format: Tubes

Shipping condition: Dry ice

Product Components:

- Glycerol-Free T4 UvsX Recombinase
- 10x UvsX Recombinase Reaction Buffer
- 20 mM Tris-acetate pH 7.8, 100 mM Potassium acetate 10 mM Magnesium acetate, 1 mM DTT

Storage:

-20°C

Additional Buffer Information:

Suggested Storage/ Dilution Buffer:

50 mM Tris-HCl, 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, pH 7.5 @ 25°C




Fig. 1:
Lane 1, Protein marker
Lane 2, UvsX Recombinase

Quality Control:

UvsX recombinase is free from detectable nuclease activities .

Notes:

Glycerol acts as a cryoprotectant and protein stabilizer when added to the storage buffer of proteins and enzymes. However, in certain circumstances, it is preferred to omit glycerol from the buffer. This includes instances where the presence of glycerol may interfere, such as lyophilization, high-throughput instruments with sensitive fluidics or primary cell cultures. For our glycerol free products, it is recommended to use immediately upon thaw as without glycerol there will be significant activity loss with freeze/thaw cycles.

Enzyme Concentration:

IG uses orthogonal, 3-part approaches to determine the enzyme concentration to provide you with consistent and reliable enzymes for your needs. The quantity of a protein sample is assessed using densitometry with polyacrylamide gel electrophoresis (PAGE), UV absorbance spectra of native protein, and using a protein standard assay such as bicinchoninic acid assay (BCA) using bovine serum albumin (BSA) as a standard (Figure 1).

Why does IG use all three approaches?

1. Each method above has limitations. The limitations include experimental noise, accuracy, and susceptibility to buffer and/or enzyme conditions.
2. Each enzyme has unique physical properties that make a single approach to analyzing proteins a challenge. Each enzyme has a different protein sequence, different requirements to be stable in solution, and different requirements to retain its maximal activity. These differences can interfere with or convolute results, especially when compared to other enzymes. When used together, however, each method provides the scientist with independent measures of both enzyme and buffer purity and quality.




Figure 1: Enzyme quantitation methods used by IG.

- A) SDS-polyacrylamide gel electrophoresis. Ladder in 1st lane, 2 μ g BSA (~67 kDa) as a standard in lanes 1-3, and IG enzymes (~40 kDa) in lanes 4-6. The yellow boxes are the areas evaluated by densitometry. The integrated band intensities of IG enzymes are compared with integrated band intensities from BSA to assay concentration.
- B) UV spectrum of a clean IG enzyme with protein peaks at 230 nm and at 280 nm. An extinction coefficient at 280 nm is typically used to quantify protein using these spectra with buffer subtraction at 330 nm.
- C) BCA standard curve for BSA. The curve is used to calculate an IG enzyme concentration using BSA as the standard.

Related Products:

- RPA Kit (Cat.#3526 & #3530)
- T4 gp32 Protein (Cat.# 3515)
- T4 UvsY Protein (Cat.# 3572)
- Bsu DNA Polymerase (Cat.# 3585)
- *Sau* DNA Polymerase (Cat.# 3595)
- Exonuclease III (Cat.# 3415)
- Exonuclease IV (Nfo) (Cat.# 3425)

Ordering Information:

- Order online within the USA. Place orders on www.intactgenomics.com using our secure Shopping Cart.
- Order by email, phone, or fax.
Email: sales@intactgenomics.com
Phone: (314) 942-3655 | Toll-free : 855-835-7172 | Fax: (314) 942-3656
- Order via our distributors.

References:

1. Cromie GA, Connelly JC, Leach DR (2001) Recombination at double-strand breaks and DNA ends: conserved mechanisms from phage to humans. *Mol Cell* 8: 1163–1174
2. Michel B, Grompone G, Flores MJ, Bidnenko V (2004) Multiple pathways process stalled replication forks. *Proc Natl Acad Sci U S A* 101: 12783–12788
3. Liu J, Ehmsen KT, Heyer WD, Morrical SW (2011) Presynaptic filament dynamics in homologous recombination and DNA repair. *Crit Rev Biochem Mol Biol* 46: 240–270

Intact Genomics owns the following registered trademarks granted by the United States Patent and Trademark Office (USPTO): Intact Genomics®, IG®, ig®, igTherapeutics®, FastAmp®, i7®, DirectPlate®.

All technology protocols discussed within this manual are assumed proprietary to Intact Genomics. This Product may be covered by pending or issued patents or may have certain limitations. Please contact us for more information. Purchase of this material conveys to buyer the non-transferable right to use the material purchased in research conducted by buyer, whether for teaching, non-commercial or commercial research purposes. Buyer may not sell or otherwise transfer these materials, its components, or unmodified descendants to a third party.

Product Use Limitation and Disclaimers

This product is for research purposes only. It is not intended for therapeutic or diagnostic purposes in humans or animals. This product contains chemicals which may be harmful if misused or direct human contact is made.

Intact Genomics is dedicated to practicing and maintaining science and technology ethics. Buyer agrees to use the purchased materials in full compliance with applicable law and regulations.

Technical Support & Customer Services

Intact Genomics (IG®) is dedicated to customer satisfaction regarding the use of our products for your research needs. Each new lot of our products is thoroughly tested to ensure it meets high quality standards and provides excellent results. We appreciate your business and your feedback regarding the performance of our products in your applications. Please follow the instructions carefully and contact us if additional assistance is needed.

Our hours are Monday - Friday, 8AM to 5PM, U.S. Central Standard Time.

Intact Genomics, Inc.

11840 Westline Industrial Drive, Suite 120,
St. Louis, MO. 63146, USA

Phone: (314) 942-3655 | **Toll-free :** 855-835-7172 | **Fax:** (314) 942-3656

Email: sales@intactgenomics.com | ig@intactgenomics.com

Website: www.intactgenomics.com

© 2024 Intact Genomics, Inc
All Rights Reserved

