



Taq DNA Polymerase 2x Premix with Dye

Manual

Catalog #	3249	3250
Package Size	500 reactions	1,000 reactions
Volume	12.5 ml	25 ml



Important!

-20°C Storage Required

- * Immediately inspect packages
- * Freeze upon receipt



visit us online for more
products & custom services

Intact Genomics, Inc.

Table of Contents

Product Description.....	3
Applications.....	3
Product Components.....	3
Storage.....	3
1x Premix Composition.....	3
Protocol.....	4
Related Products.....	5
Ordering Information.....	5
References.....	5
Technical Support.....	6

Description:

Intact Genomics (ig®) *Taq* DNA Polymerase 2x Premix with Dye is a ready to use premix which contains *Taq* DNA polymerase, dNTPs, MgCl₂ and stabilizers with optimized reaction buffer. It has been optimized for routine PCR applications. *Taq* is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity (1, 2) and a 5' flap endonuclease activity (3, 4). This product is supplied with the unique Intact Genomics 5x Magic Enhancer that enables efficient amplification of GC rich templates up to 84%.

Applications:

- Routine PCR and RT-PCR
- Primer extension
- Colony PCR
- Genotyping
- Efficient for amplifying high GC content template DNA with Magic Enhancer.

Product Components:

- *Taq* DNA Polymerase 2x Premix with Dye
- 5x Magic Enhancer

Storage:

-20 °C

1x Premix Composition:

- 10 mM Tris-HCl pH 9.0
- 50 mM KCl
- 1.5 mM MgCl₂
- 0.2 mM dNTPs
- 5% Glycerol
- 0.08% Igepal CA 630
- 0.05% Tween-20
- 50 Units/ml *Taq* Polymerase

Transformation Protocol:

1. Thaw primer solutions, 5x Magic Enhancer (if required) and mix thoroughly before use.
2. Prepare a reaction mix according to the following table:

PCR Reaction Set Up:	
Template	1-50 ng
Forward Primer (3.2 μ M)	1.0 μ l
Reverse Primer (3.2 μ M)	1.0 μ l
5x Magic Enhancer (optional)	4.0 μ l
Taq DNA Polymerase 2x Premix with Dye	10.0 μ l
H ₂ O up to	20.0 μ l

3. Mix the reaction mixture thoroughly.
4. Add template DNA to the individual PCR tubes containing the reaction mixture.
5. Program the thermal cycler according to the manufacturer's instructions. A typical PCR cycling program is outlined in the following table:

PCR Cycling Conditions			
Steps	Temp.	Time	Cycles
Initial Denaturation	95 °C	3 min	1
Denaturation	95 °C	30 sec	25-35
Annealing	55-60 °C	30 sec	
Extension	72 °C	1 min/kb	
Final Extension	72 °C	5-10 min	1
Hold	4 °C	∞	

6. Place the PCR tubes in the thermal cycler and start the cycling program.
7. Analyze 5 μ l of PCR products by agarose gel electrophoresis.

Related Products:

- Hot Start *Taq* DNA Polymerase (Cat.# 3293)
- *Taq* DNA Polymerase (Cat.# 3243)
- i7® High-Fidelity DNA Polymerase 2X Master Mix (Cat.# 3257, 3259)
- i7® Hot Start High-Fidelity DNA Polymerase 2X Master Mix (Cat.# 3284, 3286)

Ordering Information:

- Order online within the USA. Place orders on **www.intactgenomics.com** using our secure Shopping Cart.
- Order by email, phone, or fax.
Email: **sales@intactgenomics.com**
Phone: (314) 942-3655 | Toll-free : 855-835-7172 | Fax: (314) 942-3656
- Order via our distributors.

References:

1. Chien, A., Edgar, D.B. and Trela, J.M. (1976). *J. Bact.* 127, 1550-1557.
2. Lawyer, F.C. et al. (1993). *PCR Methods and Appl.* 2, 275-287.
3. Longley, M.J., Bennett, S.E. and Mosbaugh D.W. (1990). *Nucleic Acids Res.* 18, 7317-7322.
4. Lyamichev, V., Brow, M.A. and Dahlberg, J.E. (1993). *Science.* 260, 778-783.

Intact Genomics owns the following registered trademarks granted by the United States Patent and Trademark Office (USPTO): Intact Genomics®, IG®, ig®, igTherapeutics®, FastAmp®, i7®, DirectPlate®.

All technology protocols discussed within this manual are assumed proprietary to Intact Genomics. This Product may be covered by pending or issued patents or may have certain limitations. Please contact us for more information. Purchase of this material conveys to buyer the non-transferable right to use the material purchased in research conducted by buyer, whether for teaching, non-commercial or commercial research purposes. Buyer may not sell or otherwise transfer these materials, its components, or unmodified descendants to a third party.

Product Use Limitation and Disclaimers

This product is for research purposes only. It is not intended for therapeutic or diagnostic purposes in humans or animals. This product contains chemicals which may be harmful if misused or direct human contact is made.

Intact Genomics is dedicated to practicing and maintaining science and technology ethics. Buyer agrees to use the purchased materials in full compliance with applicable law and regulations.

Technical Support & Customer Services

Intact Genomics (IG®) is dedicated to customer satisfaction regarding the use of our products for your research needs. Each new lot of our products is thoroughly tested to ensure it meets high quality standards and provides excellent results. We appreciate your business and your feedback regarding the performance of our products in your applications. Please follow the instructions carefully and contact us if additional assistance is needed.

Our hours are Monday - Friday, 8AM to 5PM, U.S. Central Standard Time.

Intact Genomics, Inc.

11840 Westline Industrial Drive, Suite 120,
St. Louis, MO. 63146, USA

Phone: (314) 942-3655 | **Toll-free :** 855-835-7172 | **Fax:** (314) 942-3656

Email: sales@intactgenomics.com | ig@intactgenomics.com

Website: www.intactgenomics.com

