



igScript™ Probe-Based RT-qPCR Kit

Manual

| | | | |
|---------------------|--------------|---------------|---------------|
| Catalog # | 4243 | 4245 | 4247 |
| Package Size | 50 reactions | 100 reactions | 500 reactions |
| Volume | 20µl | | |



Important!

-20°C Storage Required

- * Immediately inspect packages
- * Freeze upon receipt



visit us online for more
products & custom services

Intact Genomics, Inc.

Table of Contents

| | |
|------------------------------------|------------|
| Product Description..... | 3 |
| Components and Storage..... | 3 |
| Applications..... | 4 |
| Benefits..... | 4 |
| Protocol..... | 4,5 |
| Related Products..... | 6 |
| Ordering Information..... | 6 |
| Technical Support..... | 7 |

Description:

igScript™ Probe-Based RT-qPCR Kit combines two powerful mixtures: i). igScript™ First Strand cDNA Synthesis Kit and ii) igScript™ Probe-Based qPCR 2x Master Mix with standard buffer providing improved PCR efficiency, wider dynamic range, superior sensitivity and specificity. The two mixtures require minimal handling during reaction setup and offer consistent and robust RT-qPCR reactions.

igScript™ First Strand cDNA Synthesis Kit includes 5x igScript™ Master Mix which contains igScript™ Reverse Transcriptase, RNase inhibitor, dNTPs, an optimized buffer, MgCl₂ and protein stabilizers. igScript™ Reverse Transcriptase is a recombinant MMLV reverse transcriptase with reduced RNase H activity and increased thermostability. The kit also provides two optimized primers and nuclease-free water. An anchored Oligo-dT primer [d(T)₂₃VN] forces the primer to anneal to the beginning of the polyA tail and the random hexamer primer mix provides random and consistent priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs. The kit is highly efficient at producing full-length cDNA from long RNA templates at temperatures between 42-50 °C.

igScript™ Probe-Based qPCR 2x Master Mix is a ready-to-use cocktail containing all components except primers and template, for the amplification and detection of DNA. This 2x master mix includes proprietary buffer, *Taq* DNA polymerase, low ROX reference dye, MgCl₂, dNTPs, stabilizers and enhancers.

Product Components:

- 5x igScript™ Master Mix
- Oligo d(t)₂₃VN primer (50 μM)
- Random Hexamer primer mix (60 μM)
- igScript™ Probe-Based qPCR 2x Master Mix
- Nuclease-free water

Storage Temperature:

-20°C

Applications:

- Gene expression data validation
- Multiplexing
- Mutation detection
- Pathogen and viral detection
- Genetically modified organisms (GMO) characterization and genetic profiling

Benefits:

- Enhanced efficiency, specificity, and sensitivity
- Compatible with all real-time PCR instruments
- Superior gene expression results under various cycling conditions
- Robust and active for cDNA synthesis at temperatures up to 50°C.

Protocol:

(A). First strand cDNA synthesis:

1. In sterile micro-centrifuge tube, add the following components on ice:

| Component | Volume |
|--|---------------|
| Total RNA | Up to 1.0 µg |
| 5x igScript™ Master Mix | 4.0 µl |
| Primer: d(T) ₂₃ VN (50 µM) and/or random primer mix (60 µM) or Gene specific primer (10 µM) | 2.0 µl |
| Nuclease-free water | Up to 20.0 µl |

2. If using random hexamers, incubate the reaction mixture at 25°C for 10 minutes, then proceed to step 3.
3. Incubate the reaction mixture at temperatures between 42°C to 50°C for 30-60 minutes.
4. Inactivate the reaction by incubating at 65°C for 20 minutes.
5. Store products at -20°C or proceed to next step.

(B). PCR Amplification:

1. Place all kit components and cDNA samples on ice.
2. Mix and then centrifuge briefly to collect contents at the bottom of the tube.
3. Prepare a master mix for each reaction and control reaction plus 10% extra to allow for pipetting error according to the following table:

| PCR Reaction Set Up: | |
|--|------------|
| Diluted cDNA | 1.0-5.0 µl |
| Forward primer (5 µM) | 1.0 µl |
| Reverse primer (5 µM) | 1.0 µl |
| Probe (5 µM) | 0.5 µl |
| igScript™ Probe-Based qPCR 2x Master Mix | 10.0 µl |
| H ₂ O up to | 20 µl |

4. Mix the reaction mixture thoroughly.
5. Program the thermal cycler according to the manufacturer's instructions.
6. A typical PCR cycling program is outlined in the following table:

| PCR Cycling Conditions | | | |
|--------------------------|------------------------------------|--------|--------|
| Steps | Temperature | Time | Cycles |
| Initial denaturation | 95°C | 3 min | 1 |
| Denaturation | 95°C | 10 sec | 40 |
| Annealing/ Extension* | 55-60°C | 30 sec | |
| Melting curve analysis | According to instrument guidelines | | |

7. Place the PCR tubes in the thermal cycler and start the cycling program.
8. Analyze the data according to the manufacturer protocol.

Note: For 3 step cycling protocols, anneal at optimal annealing temperature for 30 sec followed by the minimum time required for data acquisition at 72°C according to instrument guidelines.

Related Products:

- igScript™ One Step RT-PCR Kit (Cat.# 4211)
- igScript™ One Step RT-qPCR Kit (Cat.# 4214)
- igScript™ Probe-Based qPCR 2x master mix (Cat.# 4233, 4235, 4237)
- igScript™ Probe-Based One Step RT-qPCR Kit (Cat.#4243, 4245, 4247)
- igScript™ First Strand cDNA Synthesis Kit (Cat.# 4312)
- igScript™ Reverse Transcriptase (Cat.# 3344)
- ig® SYBR Green qPCR 2x Master Mix (Cat.# 3354)

Ordering Information:

- Order online within the USA. Place orders on **www.intactgenomics.com** using our secure Shopping Cart.
- Order by email, phone, or fax.
Email: **sales@intactgenomics.com**
Phone: (314) 942-3655 | Toll-free : 855-835-7172 | Fax: (314) 942-3656
- Order via our distributors.

Intact Genomics owns the following registered trademarks granted by the United States Patent and Trademark Office (USPTO): Intact Genomics®, IG®, ig®, igTherapeutics®, FastAmp®, i7®, DirectPlate®.

All technology protocols discussed within this manual are assumed proprietary to Intact Genomics. This Product may be covered by pending or issued patents or may have certain limitations. Please contact us for more information. Purchase of this material conveys to buyer the non-transferable right to use the material purchased in research conducted by buyer, whether for teaching, non-commercial or commercial research purposes. Buyer may not sell or otherwise transfer these materials, its components, or unmodified descendants to a third party.

Product Use Limitation and Disclaimers

This product is for research purposes only. It is not intended for therapeutic or diagnostic purposes in humans or animals. This product contains chemicals which may be harmful if misused or direct human contact is made.

Intact Genomics is dedicated to practicing and maintaining science and technology ethics. Buyer agrees to use the purchased materials in full compliance with applicable law and regulations.

Technical Support & Customer Services

Intact Genomics (IG®) is dedicated to customer satisfaction regarding the use of our products for your research needs. Each new lot of our products is thoroughly tested to ensure it meets high quality standards and provides excellent results. We appreciate your business and your feedback regarding the performance of our products in your applications. Please follow the instructions carefully and contact us if additional assistance is needed.

Our hours are Monday - Friday, 8AM to 5PM, U.S. Central Standard Time.

Intact Genomics, Inc.

11840 Westline Industrial Drive, Suite 120,
St. Louis, MO. 63146, USA

Phone: (314) 942-3655 | **Toll-free :** 855-835-7172 | **Fax:** (314) 942-3656

Email: sales@intactgenomics.com | ig@intactgenomics.com

Website: www.intactgenomics.com

© 2024 Intact Genomics, Inc
All Rights Reserved

